SHANGHAI VANA BIOTECH CO., LTD.



http://www.vanabio.com

Add:2-C03 Room,2 Build, No. 1588 of Huhang Road, Shanghai City, China.

SILIT-3300

TECHNICAL DATA SHEET

Introduction

SILIT-3300 is one block silicone softener, the product can be used as various textiles finishing agent (such as cotton and its blends, rayon, viscose fiber, synthetic fiber, silk, wool, etc). Especially suited for cotton and blends fabrics. **SILIT-3300** have super soft and good smooth handfeeling.

Characteristics

Appearance: Transparent yellow liquid pH value: 4~6(1% Aqueous solution)

Solid content: 80%

Ionic Character: Light cationic

Diluent: Water

Compatibility: Be mixed easily with cationic and nonionic auxiliaries

Features of products

- Provides super softness and good smoothness.
- Very low yellowing.
- Excellent product stability.
- Provides dilutions that have excellent compatibility with conventional textile auxiliaries.

Applications

Only one thing need be attention. In fact **SILIT-3300** is oil, it need make chemical emulsion inversion around 30% solid content by carefully stirring.

So factory must seriously stir it before using, pls strictly dilute it with the following method.

- 375kgs SILIT-3300, first add 125kgs water, keep stirring 20-30 minutes, until the emulsion is homogenous and transparent.
- 2) continue to add 250kgs water, keep stirring 10-20 minutes until the emulsion is homogenous and transparent.
- 3) continue to add 250kgs water, keep stirring 10-20 minutes until the emulsion is homogenous and transparent.

So now it is 30% solid content emulsion and stable enough

1.Exhaustion process:



SHANGHAI VANA BIOTECH CO., LTD.

http://www.vanabio.com

Add:2-C03 Room,2 Build, No. 1588 of Huhang Road, Shanghai City, China.

SILIT-3300 (30%,after 1~3%o.w.f.

dilution)

Usage :40 °C ~50 °C ×15~30min

2.Padding process:

SILIT-3300(30%,after 10~30g/L

dilution)

Usage :once-dip-once-nip

Package

SILIT-3300 is available in 200kg plastic drums.

Storage and shelf-life

When stored in its original packaging at a temperature of between -20°C and +50°C, **SILIT-3300** may be stored for up to 12 months from its date of manufacture (expiry date). Comply with the storage instructions and expiry date marked on the packaging. Past this date, **Shanghai Vana Biotech** no longer guarantees that the product meets the sales specifications.